

5-Methylcytosine & 5-Hydroxymethylcytosine DNA Standard Set



Cat. No. D5405

Storage: -20 °C


Product Information

Applications:

- Control DNA for 5-methylcytosine and 5-hydroxymethylcytosine quantitation applications (i.e. - mass spectrometry, HPLC, TLC, etc.)¹
- Substrate for studies involving 5-hydroxymethylcytosine interacting proteins

Description:

The 5-Methylcytosine & 5-Hydroxymethylcytosine DNA Standard Set is a set of three DNA standards that are linear dsDNA, 897 bp, and have the same sequence (see sequence below & figure). The only difference is that each contains either 100% unmodified cytosines, 5-methylcytosines, or 5-hydroxymethylcytosines (see figure). Since the sequence and extent of cytosine modification is known, this DNA standard set is ideal for use in calibration of various applications intended for quantitation of cytosine modifications.



5-Methylcytosine & 5-Hydroxymethylcytosine DNA Standard Set
All three standards are 897 bp with all cytosines either unmodified (lane 1, C), 5-methylcytosine (lane 2, ^{5m}C), or 5-hydroxymethylcytosine (lane 3, ^{5hm}C). 1kb DNA ladder (Zymo Research).

Product Contents:

	Cat. No.	DNA	Storage
Cytosine DNA Standard	D5405-1	2 µg	-20°C
5-Methylcytosine DNA Standard	D5405-2	2 µg	-20°C
5-Hydroxymethylcytosine DNA Standard	D5405-3	2 µg	-20°C

Storage: Store reagents at -20 °C for up to 12 months. Prolonged storage should be ≤ -70 °C.

DNA Concentration: [50 ng/µl] in storage buffer (10mM Tris-HCl; 0.1mM EDTA; pH 8.0). All standards are supplied in 40 µl volume.

Sequence Information

```

5'... 1 AGAATTGGTT AATTGGTTGT AACACTGGCA GAGCATTACG CTGACTTGAC
51 GGGACGGCGG CTTTGTGTA TAAATCGAAC TTTTGTGAG TTGAAGGATC
101 AGATCACGCA TCTTCCCGAC AACCGAGACC GTTCCGTGGC AAAGCAAAAAG
151 TTCAAAAATCA CCAACTGGTC CACCTACAAC AAAGCTCTCA TCAACCGTGG
201 CTCCTCACT TTCTGGCTGG ATGATGGGGC GATTACGGCC TGGTATGAGT
251 CAGCAACACC TTCTTCACGA GGCAGACCTC AGCGCTCAAA GATGCAGGGG
301 TAAAAGCTAA CCGCATCTTT ACCGACAAGG CATCCGGCAG TTCAACAGAT
351 CGGGAAGGGC TGGATTGGCT GAGGATGAAG GTGGAGGAAG GTGATGTCAT
401 TCTGGTGAAG AAGCTCGACC GTCTTGGCCG CGACACCGCC GACATGATCC
451 AACTGATAAA AGAGTTTGAT GCTCAGGGTG TAGCGGTTTC GTTTATTGAC
501 GACGGGATCA GTACCGACGG TGATATGGGG CAAATGGTGG TCACCATCCT
551 GTCGGCTGTG GCACAGGCTG AACGCCGGAG GATCCTAGAG CGCACGGAATG
601 AGGGCCGACA GGAAGCAAAG CTGAAAGGAA TCAAATTTGG CCGCAGGCGT
    
```

```

651 ACCGTGGACA GGAACGTCGT GCTGACGCTT CATCAGAAGG GCACTGGTGC
701 AACGGAAATT GCTCATCAGC TCAGTATTGC CCGCTCCACG GTTTATAAAA
751 TTCTTGAAGA CGAAAGGGCC TCGTGATACG CCTATTTTTA TAGGTTAATG
801 TCATGATAAT AATGGTTTCT TAGACGTACG GTGGCACTTT TCGGGGAAAT
851 GTGCGCGGAA CCCTATTTG TTTATTTTTC TAAATACATT CAAATAT...3'
    
```

Note: Sequence is same for all three standards. Depending on the DNA, all C's are unmodified cytosine, 5-methylcytosine, or 5-hydroxymethylcytosine. Sequence can be downloaded on the product page of our website.

Sequence Statistics

Nucleotide	Count	Frequency (%)
Adenine (A)	237	26.4
Cytosine (C)	201	22.4
Guanine (G)	246	27.4
Thymine (T)	213	23.7
C + G	447	49.8
A + T	450	50.2

References

- Ito, S. *et al.* (2010) "Role of Tet proteins in 5mC to 5hmC conversion, ES-cell self-renewal and inner cell mass specification". *Nature* 466, 1129-1133.

Related Epigenetics Products:

Product Name	Size	Cat. No.
5-HYDROXYMETHYL CYTOSINE		
Quest 5-hmC Detection Kit™	25 Preps. 50 Preps.	D5410 D5411
Quest 5-hmC Detection Kit™ - Lite	25 Preps. 50 Preps.	D5415 D5416
Human Matched DNA Set	2 x 5 µg	D5018
Mouse 5hmC & 5mC DNA Set	4 x 5 µg	D5019
5-hmC Glucosyltransferase	100 units 200 units	E2026 E2027
5-Hydroxymethyl dCTP [100mM]	10 µmol	D1045
5-Methyl dCTP [10mM]	1 µmol	D1035
5-Methylcytosine & 5-Hydroxymethylcytosine DNA Standard Set	1 set	D5405
BISULFITE TREATMENT OF DNA		
EZ DNA Methylation-Direct™ Kit	50 rxns.	D5020
	200 rxns.	D5021
	2 x 96 rxns.	D5022
	2 x 96 rxns.	D5023
EZ DNA Methylation-Startup™ Kit	50 rxns.	D5024
EZ Bisulfite DNA Clean-up Kit™	50 preps.	D5025
	200 preps.	D5026
	2 x 96 preps.	D5027
	2 x 96 preps.	D5028
METHYLATED/NON-METHYLATED DNA STANDARDS		
Universal Methylated DNA Standard	1 set	D5010
Universal Methylated Human DNA Standard	1 set	D5011
Universal Methylated Mouse DNA Standard	1 set	D5012
Human Methylated and Non-methylated DNA Set	1 set	D5014
Human HCT116 DKO Non-methylated DNA Standard	5 µg	D5014-1
Human HCT116 DKO Methylated DNA Standard	5 µg	D5014-2
Bisulfite Converted Universal Methylated Human DNA Standard	1 set	D5015
<i>E. coli</i> Non-methylated Genomic DNA	5 µg	D5016

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com

5-Methylcytosine & 5-Hydroxymethylcytosine DNA Standard Set



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Cat. No. D5405

Storage: -20 °C

Product Information

AMPLIFICATION OF BISULFITE CONVERTED DNA		
Zymo Taq™ DNA Polymerase	50 rxns.	E2001
	200 rxns.	E2002
Zymo Taq™ PreMix (2X concentrated)	50 rxns.	E2003
	200 rxns.	E2004
ANTIBODIES & IMMUNOPRECIPITATION		
Methylated-DNA IP Kit	10 preps.	D5101
ChIP DNA Clean & Concentrator™	50 preps.	D5201
	50 preps.	D5205
Anti-5-Methylcytosine Monoclonal Antibody (clone 10G4)	50 µg	A3001-50
	200 µg	A3001-200
METHYLTRANSFERASES		
CpG Methylase (M.SssI)	200 U	E2010
	400 U	E2011
GpC Methylase (M.CviPI)	200 U	E2014
	1000 U	E2015
DNA FRAGMENTATION		
DNA Degradase™	500 U	E2016
	2000 U	E2017
DNA Degradase Plus™	250 U	E2020
	1000 U	E2021
DNA Shearase™	50 U	E2018-50
	200 U	E2018-200
	50 U & DCC™	E2019-50
	200 U & DCC™	E2019-200
NUCLEOSOME MAPPING		
EZ Nucleosomal DNA Prep Kit	20 preps	D5220

Version 1.0.2

Trademarks and Disclaimers:

™ Trademarks of Zymo Research Corporation.

This product is for research use only and should only be used by trained professionals. It is not intended for use in diagnostic procedures. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide maximal performance and reliability.

For Technical Assistance, please contact 1-888-882-9682 or E-mail tech@zymoresearch.com. Satisfaction of all Zymo Research products is guaranteed. If you are not satisfied with this product please call 1-888-882-9682.

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com