

Restriction  
Endonuclease



## CERTIFICATE OF ANALYSIS

Source: *Arthrobacter species 7M06*.

### Supplied in:

10 mM Tris-HCl (pH 7.5), 50 mM KCl, 0.1 mM EDTA,  
7 mM 2-mercaptoethanol, 200 µg/ml BSA, 50%  
glycerol.

### Reaction Conditions:

1× SE-Buffer AbsI. Incubate at 37° C.

### 1X SE-Buffer AbsI (pH 9.0 @ 25° C):

10 mM Tris-HCl    50 mM KCl  
10 mM MgCl<sub>2</sub>    1 mM DTT

### Heat Inactivation:

Enzyme is inactivated by incubation at 65 °C for 20  
minutes.

Unit Definition: One unit is defined as the amount of  
enzyme required to digest 1 µg of pUC19SE/Dril DNA  
in 1 hour at 37° C in a total reaction volume of 50 µl.

### Quality Control Assays

Ligation: After 2-fold overdigestion with Abs I, ~90%  
of the DNA fragments can be ligated with T4 DNA Ligase  
and recut.

16-Hour Incubation: A 50 µl reaction containing 1 µg  
of DNA and 2 Units of enzyme incubated for 16 hours  
resulted in the same pattern of DNA bands as a reaction  
incubated for 1 hour.

A long incubation time may result in star activity.

Oligonucleotide Assay: No detectable degradation of  
a single-stranded and double-stranded oligonucleotide  
was observed after incubation with 1 units of restriction  
endonuclease for 3 hours.

### Enzyme Properties:

When using a buffer other than the optimal (Supplied)  
SE-Buffer, it may be necessary to add more enzymes  
to achieve complete digestion.

### Reagents Supplied with Enzyme:

10X SE Buffer AbsI.

## Abs I

Recognition  
Sequence:

CC↓TCGAGG  
GGAGCT↑CC

S

E535

50 units

1,000 u/ml

Lot:

Exp:

Store at -20°C

SE-Buffers	B	G	O	W	Y	ROSE
%Activity	75-100	10-25	0	50-75	0-10	50

37°C

65°C

AbsI

pUC19SE/  
Dril

For more details  
scan the code



Ph/F+7(383)333-6853  
info@sibenzyme.com  
www.sibenzyme.com